

PE/Cyanine5 anti-human CD181 (CXCR1)

Mouse IgG2b, kappa
8F1/CXCR1



BioLegend
From Revvity

100 tests

Artikel-Nr.: B320610 | BioLegend | Hersteller-Nr.: 320610

244,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 8F1/CXCR1
Format: PE_CYANINE5
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human CD181 (CXCR1)

Clone: 8F1/CXCR1

Isotype: Mouse IgG2b, kappa

Other_Names: CD181, CXCR1, CDw128, CD128a, IL-8RA, Type I IL-8 receptor

Size: 100 tests

Reactivity: Human

Format: PE_CYANINE5

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3003>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>