

Purified anti-human MICA/MICB

Mouse IgG2a, kappa
6D4



BioLegend
From Revvity

100 µg
Artikel-Nr.: B320902 | **BioLegend** | Hersteller-Nr.: 320902

95,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 6D4
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human MICA/MICB

Clone: 6D4

Isotype: Mouse IgG2a, kappa

Other_Names: PERB11, MIC, MHC class I chain-related protein A, MICA, PERB11.1, MHC class I chain-related protein B, MICB, PERB11.2

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3062>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>