

Purified anti-human CD206 (MMR)

Mouse IgG1, kappa
15-2



BioLegend
From Revvity

25 µg
Artikel-Nr.: B321101 | **BioLegend** | Hersteller-Nr.: 321101

52,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 15-2
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-human CD206 (MMR)

Clone: 15-2

Isotype: Mouse IgG1, kappa

Other_Names: MMR (macrophage mannose receptor), MR (mannose receptor), CD206, MRC1

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2992>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>