

Biotin anti-human CD258 (LIGHT)

Mouse IgG2b, kappa
T5-39



BioLegend
From Revvity

100 µg

Artikel-Nr.: B318704 | **BioLegend** | Hersteller-Nr.: 318704

184,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: T5-39
Format: BIOTIN
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Biotin anti-human CD258 (LIGHT)

Clone: T5-39

Isotype: Mouse IgG2b, kappa

Other_Names: LIGHT, TL5, TNFSF14, TR2 ligand, Ligand for herpesvirus entry mediator, tumor necrosis factor receptor-like-2 ligand, tumor necrosis factor ligand superfamily member 14 (TNFSF14)

Size: 100 µg

Reactivity: Human

Format: BIOTIN

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3869>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>