

Alexa Fluor® 488 anti-mouse/rat/human FOXP3

Mouse IgG1, kappa
150D



BioLegend
From Revvity

25 tests

Artikel-Nr.: B320011 | [BioLegend](#) | Hersteller-Nr.: 320011

151,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: 150D
Format: A488
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

Alexa Fluor® 488 anti-mouse/rat/human FOXP3

Clone: 150D

Isotype: Mouse IgG1, kappa

Other_Names: Forkhead box protein P3, Scurfin, JM2, IPEX, Zinc finger protein JM2

Size: 25 tests

Reactivity: Human, Mouse, Rat

Format: A488

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2891>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>