

Purified anti-human FOXP3

Mouse IgG1, kappa
259D



BioLegend
From Revvity

25 µg
Artikel-Nr.: B320201 | **BioLegend** | Hersteller-Nr.: 320201

104,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 259D
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-human FOXP3

Clone: 259D

Isotype: Mouse IgG1, kappa

Other_Names: Forkhead box protein P3, Scurfin, JM2, IPEX, Zinc finger protein JM2

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2903>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>