

Alexa Fluor® 488 anti-human FOXP3

Mouse IgG1, kappa
259D



BioLegend
From Revvity

25 tests

Artikel-Nr.: B320211 | BioLegend | Hersteller-Nr.: 320211

132,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 259D
Format: A488
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

Alexa Fluor® 488 anti-human FOXP3

Clone: 259D

Isotype: Mouse IgG1, kappa

Other_Names: Forkhead box protein P3, Scurfin, JM2, IPEX, Zinc finger protein JM2

Size: 25 tests

Reactivity: Human

Format: A488

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2908>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>