

Purified anti-human CD268 (BAFF-R)

Mouse IgG1, kappa
11C1



BioLegend
From Revvity

100 µg
Artikel-Nr.: B316902 | [BioLegend](#) | Hersteller-Nr.: 316902

165,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 11C1
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD268 (BAFF-R)

Clone: 11C1
Isotype: Mouse IgG1, kappa
Other_Names: TNFRSF13C, BAFF-R, BAFFR, BR3, BAFF Receptor
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1720>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>