

Purified anti-human CD98

Mouse IgG1, kappa
MEM-108



BioLegend
From Revvity

100 µg

Artikel-Nr.: B315602 | **BioLegend** | Hersteller-Nr.: 315602

90,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, IP
Clone: MEM-108
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD98

Clone: MEM-108
Isotype: Mouse IgG1, kappa
Other_Names: 4F2, FRP-1, RL-388
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, IP
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2787>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>