

FITC anti-human CD98

Mouse IgG1, kappa
MEM-108



BioLegend
From Revvity

25 tests

Artikel-Nr.: B315603 | [BioLegend](#) | Hersteller-Nr.: 315603

81,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC
Clone: MEM-108
Format: FITC
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

FITC anti-human CD98

Clone: MEM-108
Isotype: Mouse IgG1, kappa
Other_Names: 4F2, FRP-1, RL-388
Size: 25 tests

Reactivity: Human
Format: FITC
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2788>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>