

Purified anti-human CD354 (TREM-1)

Mouse IgG2b, kappa
TREM-37



BioLegend
From Revvity

100 µg

Artikel-Nr.: B316102 | [BioLegend](#) | Hersteller-Nr.: 316102

184,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: TREM-37
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD354 (TREM-1)

Clone: TREM-37

Isotype: Mouse IgG2b, kappa

Other_Names: Triggering receptor expressed on myeloid cells 1, triggering receptor expressed on monocytes 1

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2826>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>