

FITC anti-human beta2-microglobulin

Mouse IgG1, kappa
2M2



BioLegend
From Revvity

100 tests

Artikel-Nr.: B316304 | [BioLegend](#) | Hersteller-Nr.: 316304

184,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: 2M2
Format: FITC
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

FITC anti-human beta2-microglobulin

Clone: 2M2
Isotype: Mouse IgG1, kappa
Other_Names: beta2M, beta2-M, beta2-microglobulin b2-M, b2M
Size: 100 tests

Reactivity: Human
Format: FITC
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3079>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>