

PE anti-human Ig light chain kappa

Mouse IgG1, kappa
MHK-49



BioLegend
From Revvity

25 tests

Artikel-Nr.: B316507 | **BioLegend** | Hersteller-Nr.: 316507

88,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: MHK-49
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human Ig light chain kappa

Clone: MHK-49
Isotype: Mouse IgG1, kappa
Other_Names: Immunoglobulin light chain kappa
Size: 25 tests

Reactivity: Human
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3298>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>