

FITC anti-human CD46

Mouse IgG1, kappa
MEM-258



BioLegend
From Revvity

100 tests

Artikel-Nr.: B315304 | **BioLegend** | Hersteller-Nr.: 315304

184,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: MEM-258
Format: FITC
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

FITC anti-human CD46

Clone: MEM-258

Isotype: Mouse IgG1, kappa

Other_Names: Complement membrane cofactor protein, measles virus receptor, MCP

Size: 100 tests

Reactivity: Human

Format: FITC

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2782>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>