

## Purified anti-human/mouse CD266 (Fn14, TWEAK R)

Mouse IgG2b, kappa  
ITEM-4



**BioLegend**  
From Revvity

100 µg  
Artikel-Nr.: B314102 | **BioLegend** | Hersteller-Nr.: 314102

**184,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FA, FC, ICC, ICFC, IHC, IP, WB  
**Clone:** ITEM-4  
**Format:** PURE  
**Isotype:** Mouse IgG2b, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-human/mouse CD266 (Fn14, TWEAK R)

Clone: ITEM-4

Isotype: Mouse IgG2b, kappa

Other\_Names: Fibroblast growth factor-inducible 14 kD protein (Fn14), TWEAK receptor, TNFRSF12A

Size: 100 µg

Reactivity: Human, Mouse

Format: PURE

App\_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2664>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>