

Biotin anti-human/mouse CD266 (Fn14, TWEAK R)

Mouse IgG2b, kappa
ITEM-4



BioLegend
From Revvity

25 µg
Artikel-Nr.: B314103 | BioLegend | Hersteller-Nr.: 314103

104,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC
Clone: ITEM-4
Format: BIOTIN
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Biotin anti-human/mouse CD266 (Fn14, TWEAK R)

Clone: ITEM-4

Isotype: Mouse IgG2b, kappa

Other_Names: Fibroblast growth factor-inducible 14 kD protein (Fn14), TWEAK receptor, TNFRSF12A

Size: 25 µg

Reactivity: Human, Mouse

Format: BIOTIN

App_Abbrev: FC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2670>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>