

## Purified anti-human IgM

Mouse IgG1, kappa  
MHM-88



**BioLegend**  
From Revvity

100 µg

Artikel-Nr.: B314502 | [BioLegend](#) | Hersteller-Nr.: 314502

**81,00 €\***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** FA, FC, ICC, ICFC, IHC, IP, WB  
**Clone:** MHM-88  
**Format:** PURE  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 µg

Purified anti-human IgM

Clone: MHM-88  
Isotype: Mouse IgG1, kappa  
Other\_Names: Immunoglobulin M  
Size: 100 µg

Reactivity: Human  
Format: PURE  
App\_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:  
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2877>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>