

PE anti-human CD354 (TREM-1)

Mouse IgG1, kappa
TREM-26



BioLegend
From Revvity

100 tests

Artikel-Nr.: B314906 | [BioLegend](#) | Hersteller-Nr.: 314906

244,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: TREM-26
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD354 (TREM-1)

Clone: TREM-26

Isotype: Mouse IgG1, kappa

Other_Names: Triggering receptor expressed on myeloid cells 1, triggering receptor expressed on monocytes 1

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2771>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>