

Purified anti-human CD283 (TLR3)

Mouse IgG2a, kappa
TLR-104



BioLegend
From Revvity

100 µg

Artikel-Nr.: B315002 | BioLegend | Hersteller-Nr.: 315002

235,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: TLR-104
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD283 (TLR3)

Clone: TLR-104

Isotype: Mouse IgG2a, kappa

Other_Names: Toll like receptor 3, TLR-3, CD283

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2773>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>