

PE anti-human/mouse ZAP-70

Mouse IgG1, kappa
1E7.2



BioLegend
From Revvity

25 tests

Artikel-Nr.: B313403 | [BioLegend](#) | Hersteller-Nr.: 313403

94,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: 1E7.2
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 tests

PE anti-human/mouse ZAP-70

Clone: 1E7.2

Isotype: Mouse IgG1, kappa

Other_Names: Zeta chain associated protein kinase, SYK related tyrosine kinase

Size: 25 tests

Reactivity: Human, Mouse

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2376>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>