

Purified anti-human/mouse/rat CD278 (ICOS)

Armenian Hamster IgG
C398.4A



BioLegend
From Revvity

500 µg
Artikel-Nr.: B313502 | **BioLegend** | Hersteller-Nr.: 313502

216,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: C398.4A
Format: PURE
Isotype: Armenian Hamster IgG
Species of Origin: Armenian Hamster
Verpackung: 500 µg

Purified anti-human/mouse/rat CD278 (ICOS)

Clone: C398.4A
Isotype: Armenian Hamster IgG
Other_Names: Inducible COStimulatory molecule, H4
Size: 500 µg

Reactivity: Human, Mouse, Rat
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Armenian Hamster

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2477>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>