

Purified anti-human/mouse CD49f

Rat IgG2a, kappa
GoH3



BioLegend
From Revvity

100 µg
Artikel-Nr.: B313602 | **BioLegend** | Hersteller-Nr.: 313602

96,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: GoH3
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 100 µg

Purified anti-human/mouse CD49f

Clone: GoH3
Isotype: Rat IgG2a, kappa
Other_Names: VLA-6 alpha chain, alpha6 integrin, integrin alpha6, ITGA6
Size: 100 µg

Reactivity: Human, Mouse, Cynomolgus, Rhesus
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2604>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>