

Biotin anti-human/mouse CD49f

Rat IgG2a, kappa
GoH3



BioLegend
From Revvity

100 µg

Artikel-Nr.: B313604 | **BioLegend** | Hersteller-Nr.: 313604

175,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: GoH3
Format: BIOTIN
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 100 µg

Biotin anti-human/mouse CD49f

Clone: GoH3

Isotype: Rat IgG2a, kappa

Other_Names: VLA-6 alpha chain, alpha6 integrin, integrin alpha6, ITGA6

Size: 100 µg

Reactivity: Human, Mouse, Cynomolgus, Rhesus

Format: BIOTIN

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2605>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>