

Purified anti-human CD218a (IL-18Ralpha)

Mouse IgG1, kappa
H44



BioLegend
From Revvity

100 µg

Artikel-Nr.: B313804 | [BioLegend](#) | Hersteller-Nr.: 313804

128,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: H44
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD218a (IL-18Ralpha)

Clone: H44

Isotype: Mouse IgG1, kappa

Other_Names: IL-18 receptor alpha chain, IL-18 receptor 1, IL-18 receptor related protein, CDw218a

Size: 100 µg

Reactivity: Human, Cynomolgus, Rhesus

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2616>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>