

PE anti-human CD218a (IL-18Ralpha)

Mouse IgG1, kappa
H44



BioLegend
From Revvity

100 tests

Artikel-Nr.: B313808 | BioLegend | Hersteller-Nr.: 313808

268,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: H44
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD218a (IL-18Ralpha)

Clone: H44

Isotype: Mouse IgG1, kappa

Other_Names: IL-18 receptor alpha chain, IL-18 receptor 1, IL-18 receptor related protein, CDw218a

Size: 100 tests

Reactivity: Human, Cynomolgus, Rhesus

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2618>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>