

Purified anti-human CD266 (Fn14, TWEAK Recepto...

Mouse IgG1, kappa
ITEM-1



BioLegend
From Revvity

100 µg

Artikel-Nr.: B314002 | BioLegend | Hersteller-Nr.: 314002

184,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: Apop, FC, IHC-P, WB
Clone: ITEM-1
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD266 (Fn14, TWEAK Receptor)

Clone: ITEM-1

Isotype: Mouse IgG1, kappa

Other_Names: Fibroblast growth factor-inducible 14 kD protein (Fn14), TNFRSF12A

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, Apop, WB, IHC-P

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2665>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>