

Purified anti-human/mouse ZAP-70

Mouse IgG1, kappa
1E7.2



BioLegend
From Revvity

100 µg
Artikel-Nr.: B313402 | **BioLegend** | Hersteller-Nr.: 313402

147,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 1E7.2
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human/mouse ZAP-70

Clone: 1E7.2

Isotype: Mouse IgG1, kappa

Other_Names: Zeta chain associated protein kinase, SYK related tyrosine kinase

Size: 100 µg

Reactivity: Human, Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2375>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>