

PE anti-human CD40

Mouse IgG1, kappa
HB14



BioLegend
From Revvity

100 tests

Artikel-Nr.: B313006 | **BioLegend** | Hersteller-Nr.: 313006

168,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: HB14
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD40

Clone: HB14
Isotype: Mouse IgG1, kappa
Other_Names: BP50, TNFRSF5
Size: 100 tests

Reactivity: Human, Cynomolgus, Rhesus
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2363>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>