

Purified anti-human CD117 (c-kit)

Mouse IgG1, kappa
104D2



BioLegend
From Revvity

25 µg
Artikel-Nr.: B313201 | **BioLegend** | Hersteller-Nr.: 313201

62,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: CyTOF®, FC, ICC, IP
Clone: 104D2
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-human CD117 (c-kit)

Clone: 104D2

Isotype: Mouse IgG1, kappa

Other_Names: Stem cell factor receptor, c-kit, mast cell growth factor receptor, steel factor receptor

Size: 25 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, CyTOF®, ICC, IP

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2372>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>