

PE anti-human CD117 (c-kit)

Mouse IgG1, kappa
104D2



BioLegend
From Revvity

100 tests

Artikel-Nr.: B313204 | BioLegend | Hersteller-Nr.: 313204

186,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 104D2
Format: PE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD117 (c-kit)

Clone: 104D2

Isotype: Mouse IgG1, kappa

Other_Names: Stem cell factor receptor, c-kit, mast cell growth factor receptor, steel factor receptor

Size: 100 tests

Reactivity: Human

Format: PE

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2373>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>