

Purified anti-human CD135 (Flt-3/Flk-2)

Mouse IgG1, kappa
BV10A4H2



BioLegend
From Revvity

100 µg

Artikel-Nr.: B313302 | **BioLegend** | Hersteller-Nr.: 313302

113,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: BV10A4H2
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD135 (Flt-3/Flk-2)

Clone: BV10A4H2
Isotype: Mouse IgG1, kappa
Other_Names: FMS-like tyrosine kinase-3, FLT-3, STK-1, Flk-2
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2357>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>