

PE/Cyanine5 anti-human HLA-A,B,C

Mouse IgG2a, kappa
W6/32



BioLegend
From Revvity

100 tests

Artikel-Nr.: B311408 | [BioLegend](#) | Hersteller-Nr.: 311408

221,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: W6/32
Format: PE_CYANINE5
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human HLA-A,B,C

Clone: W6/32
Isotype: Mouse IgG2a, kappa
Other_Names: Major Histocompatibility Class I, MHC class I
Size: 100 tests

Reactivity: Human, Cynomolgus, Rhesus
Format: PE_CYANINE5
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1873>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>