

APC anti-human HLA-A,B,C

Mouse IgG2a, kappa
W6/32



BioLegend
From Revvity

25 tests

Artikel-Nr.: B311409 | [BioLegend](#) | Hersteller-Nr.: 311409

110,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: W6/32
Format: APC
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 tests

APC anti-human HLA-A,B,C

Clone: W6/32

Isotype: Mouse IgG2a, kappa

Other_Names: Major Histocompatibility Class I, MHC class I

Size: 25 tests

Reactivity: Human, Cynomolgus, Rhesus

Format: APC

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1870>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>