

Pacific Blue anti-human HLA-A,B,C

Mouse IgG2a, kappa
W6/32



BioLegend
From Revvity

25 µg
Artikel-Nr.: B311417 | **BioLegend** | Hersteller-Nr.: 311417

100,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: W6/32
Format: PB
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Pacific Blue anti-human HLA-A,B,C

Clone: W6/32
Isotype: Mouse IgG2a, kappa
Other_Names: Major Histocompatibility Class I, MHC class I
Size: 25 µg

Reactivity: Human, Cynomolgus, Rhesus
Format: PB
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=3336>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>