

PE/Cyanine5 anti-human CD357 (GITR)

Mouse IgG1, kappa
621



BioLegend
From Revvity

100 tests

Artikel-Nr.: B311608 | **BioLegend** | Hersteller-Nr.: 311608

282,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 621
Format: PE_CYANINE5
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human CD357 (GITR)

Clone: 621

Isotype: Mouse IgG1, kappa

Other_Names: Activation-inducible TNFR, AITR, Glucocorticoid-inducible TNFR, TNFRSF18

Size: 100 tests

Reactivity: Human

Format: PE_CYANINE5

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3945>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>