

Purified anti-human CD267 (TACI)

Rat IgG2a, kappa
1A1



BioLegend
From Revvity

100 µg

Artikel-Nr.: B311902 | [BioLegend](#) | Hersteller-Nr.: 311902

165,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 1A1
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 100 µg

Purified anti-human CD267 (TACI)

Clone: 1A1

Isotype: Rat IgG2a, kappa

Other_Names: TNFRSF13B, CD267, Transmembrane Activator and CAML Interactor (TACI)

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1702>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>