

PE/Cyanine5 anti-human CD10

Mouse IgG1, kappa
HI10a



BioLegend
From Revvity

100 tests

Artikel-Nr.: B312206 | [BioLegend](#) | Hersteller-Nr.: 312206

238,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: HI10a
Format: PE_CYANINE5
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE/Cyanine5 anti-human CD10

Clone: HI10a

Isotype: Mouse IgG1, kappa

Other_Names: Common acute lymphoblastic leukemia antigen (CALLA), Enkephalinase, Neutral endopeptidase, Neprilysin

Size: 100 tests

Reactivity: Human, Cynomolgus, Rhesus

Format: PE_CYANINE5

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2216>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>