

Purified anti-human CD158b/j (KIR2DL2/L3/S2)

Mouse IgG2a, kappa
DX27



BioLegend
From Revvity

100 µg
Artikel-Nr.: B312602 | **BioLegend** | Hersteller-Nr.: 312602

90,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: DX27
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD158b/j (KIR2DL2/L3/S2)

Clone: DX27

Isotype: Mouse IgG2a, kappa

Other_Names: CD158b1 (KIR2DL2, p58.2), CD158b2 (KIR2DL3, p58.3), CD158j (KIR2DS2, p50.2), KIR-NKAT2, CD158b/CD158j

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2271>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>