

FITC anti-human CD158b/j (KIR2DL2/L3/S2)

Mouse IgG2a, kappa
DX27



BioLegend
From Revvity

25 tests

Artikel-Nr.: B312603 | **BioLegend** | Hersteller-Nr.: 312603

81,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: DX27
Format: FITC
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 tests

FITC anti-human CD158b/j (KIR2DL2/L3/S2)

Clone: DX27

Isotype: Mouse IgG2a, kappa

Other_Names: CD158b1 (KIR2DL2, p58.2), CD158b2 (KIR2DL3, p58.3), CD158j (KIR2DS2, p50.2), KIR-NKAT2, CD158b/CD158j

Size: 25 tests

Reactivity: Human

Format: FITC

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2280>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>