

## PE/Cyanine5 anti-human CD69

Mouse IgG1, kappa  
FN50



**BioLegend**  
From Revvity

25 tests

Artikel-Nr.: B310907 | [BioLegend](#) | Hersteller-Nr.: 310907

**104,00 €\***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** FN50  
**Format:** PE\_CYANINE5  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 25 tests

PE/Cyanine5 anti-human CD69

Clone: FN50

Isotype: Mouse IgG1, kappa

Other\_Names: Very Early Activation Antigen (VEA), Activation inducer molecule (AIM)

Size: 25 tests

Reactivity: Human

Format: PE\_CYANINE5

App\_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1673>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>