

APC/Cyanine7 anti-human CD69

Mouse IgG1, kappa
FN50



BioLegend
From Revvity

100 tests

Artikel-Nr.: B310914 | [BioLegend](#) | Hersteller-Nr.: 310914

253,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: FN50
Format: APC_CYANINE7
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

APC/Cyanine7 anti-human CD69

Clone: FN50

Isotype: Mouse IgG1, kappa

Other_Names: Very Early Activation Antigen (VEA), Activation inducer molecule (AIM)

Size: 100 tests

Reactivity: Human

Format: APC_CYANINE7

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1917>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>