

Alexa Fluor® 700 anti-human CD69

Mouse IgG1, kappa
FN50



BioLegend
From Revvity

25 µg

Artikel-Nr.: B310921 | [BioLegend](#) | Hersteller-Nr.: 310921

100,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: FN50
Format: A700
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Alexa Fluor® 700 anti-human CD69

Clone: FN50

Isotype: Mouse IgG1, kappa

Other_Names: Very Early Activation Antigen (VEA), Activation inducer molecule (AIM)

Size: 25 µg

Reactivity: Human

Format: A700

App_Abbrev: FC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3425>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>