

Purified anti-human CD24

Mouse IgG2a, kappa
ML5



BioLegend
From Revvity

25 µg
Artikel-Nr.: B311101 | **BioLegend** | Hersteller-Nr.: 311101

48,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: ML5
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-human CD24

Clone: ML5
Isotype: Mouse IgG2a, kappa
Other_Names: Ly-52, Heat Stable Antigen (HSA), Nectadrin, BA-1
Size: 25 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1806>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>