

Purified anti-human HLA-A,B,C

Mouse IgG2a, kappa
W6/32



BioLegend
From Revvity

100 µg

Artikel-Nr.: B311402 | BioLegend | Hersteller-Nr.: 311402

72,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: W6/32
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human HLA-A,B,C

Clone: W6/32

Isotype: Mouse IgG2a, kappa

Other_Names: Major Histocompatibility Class I, MHC class I

Size: 100 µg

Reactivity: Human, Cynomolgus, Rhesus

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1874>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>