

PE anti-human CD45RB

Mouse IgG2b, kappa
MEM-55



BioLegend
From Revvity

100 tests

Artikel-Nr.: B310204 | **BioLegend** | Hersteller-Nr.: 310204

215,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: MEM-55
Format: PE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 Tests

PE anti-human CD45RB

Clone: MEM-55
Isotype: Mouse IgG2b, kappa
Other_Names: CD45RB
Size: 100 tests

Reactivity: Human
Format: PE
App_Abbrev: FC, ICFC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1608>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>