

Purified anti-human CD193 (CCR3)

Mouse IgG2b, kappa
5E8



BioLegend
From Revvity

100 µg

Artikel-Nr.: B310702 | **BioLegend** | Hersteller-Nr.: 310702

147,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 5E8
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD193 (CCR3)

Clone: 5E8

Isotype: Mouse IgG2b, kappa

Other_Names: CC CKR3, MIP1-alpha receptor like-2, eotaxin receptor, CD193, CCR3

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1660>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>