

Purified anti-human CD154 (CD40L)

Mouse IgG1, kappa
24-31



BioLegend
From Revvity

100 µg
Artikel-Nr.: B310802 | **BioLegend** | Hersteller-Nr.: 310802

165,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 24-31
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD154 (CD40L)

Clone: 24-31
Isotype: Mouse IgG1, kappa
Other_Names: CD40L, gp39, TRAP, T-BAM, TNFSF5
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1664>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>