

## PE/Cyanine5 anti-human CD154 (CD40L)

Mouse IgG1, kappa  
24-31



**BioLegend**  
From Revvity

100 tests

Artikel-Nr.: B310808 | BioLegend | Hersteller-Nr.: 310808

**264,00 €\***

\*zzgl. MwSt. zzgl. Versandkosten

### Beschreibung

**App\_Abbrev:** FC, ICFC  
**Clone:** 24-31  
**Format:** PE\_CYANINE5  
**Isotype:** Mouse IgG1, kappa  
**Species of Origin:** Mouse  
**Verpackung:** 100 Tests

PE/Cyanine5 anti-human CD154 (CD40L)

Clone: 24-31  
Isotype: Mouse IgG1, kappa  
Other\_Names: CD40L, gp39, TRAP, T-BAM, TNFSF5  
Size: 100 tests

Reactivity: Human  
Format: PE\_CYANINE5  
App\_Abbrev: FC, ICFC  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:  
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1667>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>