

Purified anti-Phosphotyrosine

Mouse IgG2b, kappa
PY20



BioLegend
From Revvity

25 µg

Artikel-Nr.: B309301 | [BioLegend](#) | Hersteller-Nr.: 309301

44,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: PY20
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 25 µg

Purified anti-Phosphotyrosine

Clone: PY20
Isotype: Mouse IgG2b, kappa
Other_Names: p-Tyr
Size: 25 µg

Reactivity: Human, Mouse, Rat, All Species
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=1627>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>