

Purified anti-human CD275 (B7-H2, ICOSL)

Mouse IgG2b, kappa
2D3



BioLegend
From Revvity

100 µg

Artikel-Nr.: B309402 | **BioLegend** | Hersteller-Nr.: 309402

90,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 2D3
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD275 (B7-H2, ICOSL)

Clone: 2D3

Isotype: Mouse IgG2b, kappa

Other_Names: CD275, ICOSL, B7H2, ICOS-L, B7RP-1, B7-RP1

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1516>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>